大腸癌篩檢測試 Colorectal Cancer Screening Test

Q www.hicircle.com/coloclear

什麼是ColoClear?

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What is ColoClear?

ColoClear by Circle is a non-invasive, FIT-DNA screening test for detecting early signs of colorectal cancer using at-home sample collection. The test is designed to make colorectal cancer screening **simple** and **convenient**, so you can be proactive and in control of your health.

It is the only test for colorectal cancer screening approved by the National Medical Products Administration (NMPA).¹

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ColoClear is clinically validated in Asian population studies.

您的ColoClear測試結果代表甚麼? What do your ColoClear results mean?

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A positive result indicates that a signal for colorectal cancer lesions have been detected in the stool sample. You will likely need to plan a follow-up colonoscopy with your healthcare provider.

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A negative result indicates that based on the sample, no colorectal cancer lesions have been detected, and there is a low possibility of colorectal cancer or precancerous polyps.

您的醫護人員供者將會根據結果,為您提供下一步的建議。 Your healthcare provider will interpret your results and advise on next steps.

¹People's Republic of China. National Medical Products Administration. (2020). 體外診斷試劑產品註冊技術審評報告: CSZ2000050. https://www.cmde.org. cn/xwdt/shpbg/20201112105300335.html

早期篩檢,是您保護自己的第一步 Early screening is your first line of defense

*Survival rate here refers to the percentage of diagnosed individuals who have survived 5 years or more after their diagnosis.²

只要及早發現,高達90%的患者 都可以治療2 90% of colorectal cancer cases are treatable if detected at an early stage²

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The importance of colorectal cancer screenina

- Colorectal cancer is considered the most preventable cancer. Almost half of detected cases are found in a late stage with a low survival rate. However if it is found at an early stage, colorectal cancer has a survival rate of over 90%.²
- Cancerous polyps develop slowly over time, often showing no symptoms until a late stage.
- The risk of colorectal cancer increases significantly with age.

² United States of America. National Cancer Institute. (2022). Cancer Stat Facts: Colorectal Cancer. https://seer.cancer.gov/statfacts/html/colorect.html

掌握大腸健康的第一步 Your first step to colon health

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在進行測試前,用家應諮詢醫護專業人士意見。 Patients should consult their primary healthcare provider before conducting the test. The test is especially suitable for:

- Individuals eager to understand their colon health status and detect colorectal cancer risks early.
- Individuals above 18 years old looking for an accessible and convenient colorectal screening test.
- Individuals with family or medical history of colorectal cancer.
- Individuals with bowel symptoms, such as chronic constipation or bloating.

ColoClear有甚麼優勝之處? Why choose ColoClear for screening?

採用新一代DNA技術,多靶點糞便基因 (FIT-DNA)測試技術;技術發展成熟, 媲美美國市場領導者 Utilises next-generation DNA technology, FIT-DNA technology; a well-established technology on par with the US market leader

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Accurate results with a 99.6% negative predictive value and 96% sensitivity, clinically validated in Asian population studies

唯一獲國家藥品監督管理局 批准的大腸癌篩檢測試¹ The only approved colorectal cancer screening test by the National Medical Products Administration (NMPA)¹

非入侵性及無痛, 安在家中的採樣過程 Non-invasive and painless at-home sample collection at your convenience

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What is Stool DNA Technology?

This advanced DNA technology ensures highly accurate results by analysing multiple genetic mutations that contribute to the complex development of colorectal cancer. This way, pre-cancerous polyps and signs of colorectal cancer can be effectively detected at various stages. The Stool DNA technology is based on studies looking at DNA mutations and methylations patterns. Multiple studies have compared the genetic make up of malignant growths in colorectal cancer patients with the cells shed by healthy individuals. The results from these studies have indicated a strong correlation between certain DNA biomarkers and an increased risk for developing colorectal cancer.³

What is FIT Technology?

A traditional FIT tests for traces of blood, which can be caused by colorectal cancer and/or polyps. However, polyps may not always bleed and for this reason the sensitivity of a FIT test for colorectal cancer ranges from 73% to 82%⁴ when compared to a colonoscopy. ColoClear可檢測癌症及癌前 病變瘜肉 ColoClear detects both cancerous and pre-cancerous polyps.

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Painless at-home collection of stool sample

DA 📩

安坐家中採集樣本 Sample collection in the comfort of your home

5個工作天內獲取報告 Receive report in 5 business days 大腸鏡檢查 Colonoscopy

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入侵性大腸 內窺鏡 Invasive internal visualisation of the colon

Ð

須調整飲食2至3天 Requires adjustments in diet for 2-3 days

檢查前一晚需要 服食瀉藥 Need to take laxatives and fast night before procedure

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ColoClear is not a diagnostic test and is designed for screening purposes only. The result cannot be interpreted as a diagnosis for colorectal cancer. ColoClear can be an alternative for a colonoscopy that is intended for screening purposes. Please consult your primary healthcare provider to learn more.

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Avoid collecting the sample if you are experiencing any of the following:

- Menstrual bleeding
- Diarrhoea and/or watery stool
- Rectal bleeding (e.g. actively bleeding haemorrhoids)

ColoClear測試套裝包含甚麼? What's included in the ColoClear kit?

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本港錄得的大腸癌死亡人數 佔全部癌症的14.6% It makes up around 14.6% of all cancer deaths in Hong Kong 瘜肉有機會在無任何症狀下 演變為癌症 Polyps can grow and develop into cancer without showing symptoms

⁵Hong Kong SAR. Centre for Health Protection, Department of Health. (2022). Colorectal Cancer. https://www.chp.gov.hk/en/healthtopics/ content/25/51.htmlhealthtopics/content/25/51.html⁶Hong Kong SAR. Centre for Health Protection, Department of Health. (2020). Non-Communicable Diseases Watch – Fighting Colorectal Cancer. https://www.chp.gov.hk/files/pdf/ncd_watch_january_2020.pdf⁷Pan J. et al. (2020). Prevalence and risk factors for colorectal polyps in a Chinese population: a retrospective study. Scientific Reports. 10(1). doi: 10.1038/s41598-020-63827-6

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Where to get ColoClear?

If you'd like to get a ColoClear testing kit, please contact one of the partner clinics on our website to get your test from a healthcare provider.

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About ColoClear

ColoClear by Circle is brought in partnership between Prenetics and New Horizon Health. The technology was developed by New Horizon Health, a biotechnology pioneer in early detection cancer screening tests, aimed to promote innovation in cancer screening.

Prenetics is a global diagnostics and genetic testing company with the mission to bring health closer to people and decentralise healthcare. Prenetics has a strong scientific and medical leadership team that oversees the products developed by the company in internationally accredited certified laboratories, committed to continously developing new knowledge in different areas.

Prenetics

瀏覽以下內容,踏出保護自己的第一步 Visit our website and build your first line of defence now!

(Q www.hicircle.com/coloclear